Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK



Accredited to ISO 15189:2022

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

90 Priestly Road

The Surrey Research Park

Guildford

GU2 7AU

Contact: Dayan Wijesinghe Tel: +44 (0)1483 450 388

Fax: +44 (0)1483 565 007

E-Mail: dayan.wijesinghe@ctuk.eurofins.com

Website: https://www.eurofins.co.uk/county-pathology/

Testing performed at the above address only

DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS	Biochemistry examination activities for the purpose of clinical diagnosis	In house documented methods incorporating manufacturers' instructions as required
Serum SST *Capillary serum SST +Capillary serum (plain)	AMH CA125 CEA Cortisol DHEA-S Ferritin FSH Free PSA Free T3 Free T4 Folate HCG IgE - Total LH *+Oestradiol PSA *+Progesterone Prolactin SHBG *+Testosterone Thyroglobulin Antibody Thyroxine – Total TSH Vitamin B12 Vitamin B12 Vitamin D	Roche Cobas Modular 8000 (e801) using Electrochemiluminescence Immunoassay (ECLIA) SOP [Cobas 8000 SOP]
Saliva	Quantitative determination of salivary cortisol	

Assessment Manager: LH6 Page 1 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purpose of clinical diagnosis (cont'd)	In house documented methods incorporating manufacturers' instructions as required
Serum/Plasma *Capillary serum *Table 1.5	*Albumin *Alkaline Phosphatase *ALT AST *Bilirubin – Total C-Reactive Protein *Calcium Ceruloplasmin *Cholesterol Copper *C reactive protein (CRP) *hs-CRP *Creatinine Creatinine Kinase *GGT Glucose *HDL *Iron LD *Magnesium *Pancreatic α amylase *Phosphate *Protein – Total*Transferrin Triglyceride TIBC *Urea *Uric Acid Zinc	Roche Cobas Modular 8000 (c702) using Spectrophotometry SOP [Cobas 8000 SOP]

Assessment Manager: LH6 Page 2 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purpose of clinical diagnosis (cont'd)	In house documented methods incorporating manufacturers' instructions as required
Capillary Sample Plasma – Lithium Heparin	Albumin ALP ALT Calcium Cholesterol Creatinine CRP FT3 FT4 HDL Sodium Bilirubin Total protein Tryglyceride TSH Uric acid Urea	Roche Cobas Modular 8000 (c702) using Spectrophotometry SOP [Cobas 8000 SOP]
Venous – Lithium Heparin	Albumin ALP ALT AST Calcium Cholesterol Creatinine Creatinine Kinase CRP Ferritin FT3 FT4 GGT HDL Iron Sodium Total PSA Total protein Tryglyceride TSH Uric acid Urea	Roche Cobas Modular 8000 (c702) using Spectrophotometry SOP [Cobas 8000 SOP]

Assessment Manager: LH6 Page 3 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purpose of clinical diagnosis (cont'd)	In house documented methods incorporating manufacturers' instructions as required
		Roche Cobas Modular 8000 (c702) using Spectrophotometry SOP [Cobas 8000 SOP]
Urine	Albumin Creatinine	
Serum *Capillary serum	*Sodium Potassium NT-proBNP Lipoprotein A Free PSA	
Venous plasma/serum	Lipoprotein - associated phospholipase - A2	
	Calculated tests: Globulin LDL Cholesterol Non-HDL Cholesterol Chol/HDL ratio eGFR Corrected Calcium	
Serum	Glucose	Abbott Alinity- Spectrophotometry/Enzymatic CDx SOP-7
Serum	Albumin Pancreatic Amylase AST ALT Bilirubin Alkaline Phosphatase Calcium Magnesium Phosphate Urea Uric Acid Creatinine-Kinase Lactate-Dehydrogenase Gamma-Glutyl-Transferase Creatinine Cholesterol	

Assessment Manager: LH6 Page 4 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used	
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purpose of clinical diagnosis (cont'd)	In house documented methods incorporating manufacturers' instructions as required	
	High Density Lipoprotein Triglyceride UIBC Iron	Abbot Alinity – Calormetric CDx SOP-7	
	Carcinoembryonic antigen CA 125 Estradiol Dehydroepiandrosterone sulfate (DHEA-S) Prolactin Cortisol Leutenizing Hormone β-Human chorionic gonadotropin Follicle Stimulating Hormone Active-B12 Progesterone SHBG Folate, Testosterone Total-Prostate specific antigen, 25-OH Vitamin D assay anti-thyroid peroxidase anti-thyroglobulin Vitamin-B12 Free-Thyroxine-3 Free-Thyroxine-4 Total-Thyroxine-4 Total-Thyroxine-4 Thyroid-stimulating-hormone Free-Prostate-specific antigen NT-pro-BNP Ferritin	Abbott Alinity-Chemiluminescence immunoassay CDx SOP-7	
	Sodium Potassium	Abbott Alinity- ISE CDx SOP-7	
	CRP HS-CRP Total IGE Lipoprotein-A Ceruloplasmin Transferrin	Abbott Alinity - Turbidimetry CDx SOP-7	

Assessment Manager: LH6 Page 5 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purpose of clinical diagnosis (cont'd)	In house documented methods incorporating manufacturers' instructions as required
Faeces	Faecal haemoglobin concentration (FIT testing)	Integrated Sphere Latex Turbidimetry using HM-JACKarc II Analyser COUNT-SOP-14
Urine	Urine analysis dipstick and microscopy	Sysmex UC 3500 and UF-5000 Sysmex UF-5000 & UC-3500 SOP
Urine	Quantification of urine pH	Sysmex UC 3500
	Semi-quantification of Bilirubin, Blood, Glucose, Urine Ketones, Protein and Urobilinogen	
Urine	Quantification of white blood cells and red blood cells	Sysmex UF-5000
	Semi-quantification of epithelial cells, casts and bacteria	

Assessment Manager: LH6 Page 6 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Haematology examination activities for the purpose of clinical diagnosis	In house documented methods incorporating manufacturers' instructions as required
Venous whole blood Capillary whole blood	Full Blood Count: Basophils Eosinophils Haemoglobin Lymphocytes MCH MCHC MCV Monocytes Neutrophils HCT Platelets RBC WBC RDW MPV	Sysmex XN2000 COUNT-SOP-89
Venous whole blood Capillary whole blood	ABO and RH(D) Forward and Reverse Blood Grouping	Diamed ID card system COUNT-SOP-56
Venous whole blood Capillary whole blood	Rhesus Phenotyping	Biorad manual ID card system COUNT-SOP-56
Venous blood in EDTA	ESR	Alifax TEST 1 COUNT-SOP-64
Venous whole blood Capillary whole blood	Blood cell morphology including: White cell differential Incidental Malaria and blood borne parasites	Blood film examination using Wright Giemsa stain COUNT-SOP-72
Venous whole blood	Coagulation screen: PT APTT Fibrinogen	Stago STA Compact Max 3 coagulometer COUNT-SOP-86
Venous whole blood Capillary whole blood	HbA1c	HPLC using Tosoh G11 SOP Cdx-SOP 6

Assessment Manager: LH6 Page 7 of 8



Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Eurofins Clinical Diagnostics UK Limited

Issue No: 021 Issue date: 04 November 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Serology examination activities for the purpose of clinical diagnosis	In house documented methods incorporating manufacturers' instructions as required
Serum only	Hepatitis B Surface Antigen, Hepatitis B Surface Antibody IgG, Hepatitis B Core Total Antibody, Rubella IgG Syphillis IgG	Roche Cobas 801 COUNT-SOP-7 COBAS 8000
Venous serum/plasma	HIV-1 & HIV-2 antibodies & HIV-1 p24 antigen	Roche Cobas 801 COUNT-SOP-7 COBAS 8000 Elecsys HIV Duo
	Anti-HCV antibodies	Elecsys Anti-HCV II
END		

Assessment Manager: LH6 Page 8 of 8