


Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 <p>UKAS MEDICAL 10170</p> <p>Accredited to ISO 15189:2022</p>	<h3>Genomics England Limited</h3> <p>Issue No: 014 Issue date: 24 April 2026</p>	
	<p>Level 21 One Canada Square Canary Wharf London E14 5AB</p>	<p>Contact: Stephen Street-Howard Tel: +44 (0)808 2819 535 E-Mail: Stephen.Street-Howard@genomicsengland.co.uk Website: https://www.genomicsengland.co.uk/</p>
<p>Testing performed at the above address only</p>		

DETAIL OF ACCREDITATION

Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
Whole human genome sequencing data from an external source	<p><u>Bioinformatics</u></p> <p>Identification of Somatic and Germline regions in cancer patients which contain potentially actionable variants as follows:</p> <p>Germline Single Nucleotide Variants (SNVs) including variants with potential drug toxicity implications Indels Copy Number Variants (CNVs) larger than 2kb in size</p> <p>Somatic variants and variants of uncertain origin (tumour-normal and tumour only) Single nucleotide variants (SNVs) Indels Copy number variants (CNVs) Structural Variant (SVs) incl.</p> <ul style="list-style-type: none"> • SSX fusions • Internal tandem duplications (ITDs) • Dux4 rearrangements <p>Pan genomic biomarkers Tumour Mutational Burden</p>	<p>Documented in-house methods:</p> <p>Cancer pipeline</p> <p>Flexible scope limited to the application of the methods listed in GUI-BIO-010 and POL-BIO-003 for the agreed boundaries of software updates to the pipeline components, or the addition of new components to the pipeline which supplement existing functionality.</p>



10170
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Genomics England Limited
Issue No: 014 Issue date: 24 April 2026

Testing performed at main address only

Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>Whole human genome sequencing data from an external source (cont'd)</p>	<p><u>Bioinformatics</u> (cont'd)</p> <p>Whole human genome sequence data from an external source Identification of variants with potential diagnostic utility in genomes of individuals or families with rare diseases</p> <p>Single Nucleotide Variants (SNVs) Indels Short tandem repeats (STRs) Copy Number Variants (CNVs) larger than 2kb in size</p>	<p>Documented in-house methods:</p> <p>Rare Disease pipeline</p> <p>Flexible scope limited to the application of the methods listed in GUI-EXT-031 for the agreed boundaries of software updates to the pipeline components, or the addition of new components to the pipeline which supplement existing functionality.</p>



10170
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Genomics England Limited
Issue No: 014 Issue date: 24 April 2026

Testing performed at main address only

Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>Whole human genome sequencing data from an external source (cont'd)</p>	<p><u>Bioinformatics</u> (cont'd)</p> <p>Identification of Somatic and Germline regions in cancer patients which contain potentially actionable variants as follows:</p> <p>Germline variants (tumour normal) Single Nucleotide Variants (SNVs) including variants with potential drug toxicity implications Indels Copy Number Variants (CNVs) larger than 2kb in size</p> <p>Somatic variants and variants of uncertain origin (tumour-normal and tumour only) Single nucleotide variants (SNVs) Indels Copy number variants (CNVs) Structural Variant (SVs) incl.</p> <ul style="list-style-type: none"> • SSX fusions • Internal tandem duplications (ITDs) • Dux4 rearrangements <p>Pan genomic biomarkers Tumour Mutational Burden</p>	<p>Documented in-house methods:</p> <p>Alignment, Variant calling, QC steps and checks and identification of potentially clinical actionable variation from WGS data from cancer samples</p> <p>Germline genome alignment Dragen aligner</p> <p>Somatic genome alignment Dragen aligner</p> <p>Germline variant callers Dragen variant caller (small variants and copy number variants >2kb)</p> <p>Somatic variant callers Strelka (SNVs and indels) Somatic copy number variants (CNVs) using DRAGEN (CNVs) Manta (SVs) JuLI (SSX fusions) Pindel (ITDs) Somatic DUX4 rearrangements using Pelops (DUX4) rearrangements</p> <p>Software and platforms Nextflow Workflow configuration settings for the Genie workflow selection service (WSS) Open CGA (LIMS) Cellbase (Genomic reference and annotation database) PanelApp (Knowledgebase of gene phenotype association) Interpretation Portal Decision Support Software (DSS)</p> <p>SOPs: GUI-BIO-010 Cancer Genome Analysis Guide POL-BIO-003 Cancer Genome Bioinformatics Pipeline</p>



10170
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Genomics England Limited
Issue No: 014 Issue date: 24 April 2026

Testing performed at main address only

Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>Whole human genome sequencing data from an external source (cont'd)</p>	<p><u>Bioinformatics</u> (cont'd)</p> <p>Identification of variants (as listed below) with potential diagnostic utility in genomes of individuals or families with rare diseases.</p> <p>Single Nucleotide Variants (SNVs) Indels Short tandem repeats (STRs) Copy number variants (CNVs) larger than 2kb in size</p>	<p>Documented in-house methods:</p> <p>Alignment, Variant Calling, QC steps and checks and identification of potentially clinically actionable variation from WGS data from affected individuals with rare disease and their family members</p> <p>Genome alignment Dragen Aligner</p> <p>Variant Callers Dragen variant caller (SNVs/Indels) Dragen/ExpansionHunter (STRs) Dragen CNV (CNVs >2kb) and DragenSV (CNVs 2-10kb)</p> <p>Software and platforms Nextflow Workflow configuration settings for the Genie workflow selection service (WSS) Family selection workflow (assurance check) Genetic v reported (assurance check) Open CGA (LIMS) Exomiser (Variant prioritisation) PanelApp (Knowledgebase of gene phenotype association) Cellbase (Genomic reference and annotation database) Interpretation Portal Congenica (Decision Support Software) Clinical Variant Ark (CVA, Knowledgebase) Interpretation Browser (Decision Support Software)</p> <p>SOPs: GUI-EXT-031 Rare Disease Genome Analysis Guide – Online</p>
END		