## **Schedule of Accreditation**

issued by

## **United Kingdom Accreditation Service**

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK



27641

Accredited to ISO 15189:2022

#### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

41 Foley Street Contact: Husam Shlebak London Tel: +44 204 525 8805

W1W 7TS E-Mail: husam@clinilabs.co.uk
Website:www.clinilabs.co.uk

Testing performed at the above address only

#### Locations covered by the organisation and their relevant activities

#### **Laboratory locations:**

Location details		Activity
41 Foley Street London	Local contact: As above	Venepuncture
W1W 7TS		

#### **DETAIL OF ACCREDITATION**

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS	Biochemistry examination activities for the purposes of clinical diagnosis	Documented in house methods and manufacturer's instructions for:  DxC700 Beckman coulter Photometric method principle SOP: CL-CHEM -04
Blood- Serum	Albumin ALP ALT AST Gamma GT Total bilirubin Direct bilirubin Total Protein Globulin Sodium Potassium Chloride Bicarbonate Urea Creatinine eGFR Uric Acid Calcium	

Assessment Manager: HL Page 1 of 7



#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purposes of clinical diagnosis (cont'd)	Documented in house methods and manufacturer's instructions for DxC700 Beckman coulter cont'd
Blood- Serum	Magnesium Phosphate Iron UIBC TIBC Transferrin Cholesterol HDL Cholesterol LDL Cholesterol Apolipoprotein A1 Apolipoprotein B C3 C4 CK-NAC CKMB C-Reactive protein Glucose Haptoglobin Haemoglobin A1C IGA IGG IGM IGE Lactate Ceruloplasmin Amylase LDH Triglycerides LIH Calcium Creatinine Urea Magnesium Phosphate Glucose Amylase Urate	

Assessment Manager: HL Page 2 of 7



#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purposes of clinical diagnosis (cont'd)	Documented in house methods and manufacturer's instructions for
		Access-2 Beckman Coulter Immunossay Method SOP: CL- IMM- 05
Blood- Serum	Alpha feta protein AMH C153 CEA Cortisol C-Peptide DHEA-S Ferritin Folate FSH FT3 FT4 GI199 Insulin Luteinising hormone OV125 Prolactin Progesterone PSA PTH SHBG Oestradiol Testosterone Free Testosterone Free androgen index Thryoglobulin antibody Troponin I TSH Vitamin B12 Vitamin D ACE HGH	

Assessment Manager: HL Page 3 of 7



#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Biochemistry examination activities for the purposes of clinical diagnosis cont'd	Documented in house methods and manufacturer's instructions for:
Blood - serum	Electrophoresis (Paraprotein Typing and Quantification) Immunofixation Immunotyping	Capillarys 3/Hydrasys 2 Sebia Electrophoresis Protein zone electrophoresis – agarose gel SOP: CL-EP-09

Assessment Manager: HL Page 4 of 7



#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Haematological examination activities for the purposes of clinical diagnosis	Procedures documented in manufacturer's equipment manuals in conjunction with documented in- house procedures by the following methods:
EDTA whole blood	Haematology Full Blood Count (FBC) analysis of the following parameters: RBC HGB HCT MCV MCH MCHC RDW PLT MPV WBC Neutrophils Lymphocytes Monocytes Eosinophils Basophils Reticulocyte count	DxH900 and DXH560 Beckman Coulter Coulter principle with traditional method - electrical impedance. SOP: CL-HEAM -01

Assessment Manager: HL Page 5 of 7



#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Haematological examination activities for the purposes of clinical diagnosis cont'd	Procedures documented in manufacturer's equipment manuals in conjunction with documented in- house procedures by the following methods:
Blood	Erythrocyte Sedimentation Rate (ESR)	Alifax Quantitative Capillary Photometry and sedimentation method SOP: CL-HEAM- 02
Plasma (Citrate)	Routine Coagulation testing for: Prothrombin Time PT (INR) Activated Partial Thromboplastin Time (APTT) Fibrinogen (Clauss) DDimer (DD)	Werfen-ACL Top 350 Turbidimetric and Absorbance Measurement SOP: CL-COAG - 03
Blood	Adult Haemoglobinopathy analysis by detection of haemoglobin variants: Hb A Hb A2 Hb F Hb S Hb C Hb D Punjab Hb E	Sebia Octa3 using Capillary Electrophoresis SOP: CL-EP-09

Assessment Manager: HL Page 6 of 7



#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### **Clinilabs Limited**

Issue No: 001 Issue date: 20 November 2024

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	Blood Transfusion examinations to assist in clinical investigations	Procedures documented in manufacturer's equipment manuals in conjunction with documented in- house procedures by the following methods:
Blood Plasma	Blood Group by detection and identification of antigens:  A Rh D Positive AB Rh D Positive B Rh D Positive O Rh D Positive A Rh Negative AB Rh D Negative B RhD Negative O Rh D Negative O Rh D Negative Antibody Screen by detection of antibodies to:  Rh -C,D,E,c,e, Cw Kell -K, K, Kpa Duffy -Fya, Fyb MNSs- M, N, S, s Kidd -JKa, Jkb	Manual - Grifols Testing system and Reader Gel system SOP – CL-BT- 07
END		

Assessment Manager: HL Page 7 of 7