

Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 UKAS MEDICAL 31262 Accredited to ISO 15189: 2022	ARC Regulatory Limited Issue No: 002 Issue date: 18 September 2026 <table border="1"> <tr> <td data-bbox="395 439 842 716"> Unit 1 Concourse building 2 Catalyst Queen's Road Belfast BT3 9DT </td><td data-bbox="842 439 1493 716"> Contact: Seamus Watson E-Mail: Seamus.W@arc-regulatory.com Website: www.arc-regulatory.com </td></tr> </table> Testing performed at the above address only	Unit 1 Concourse building 2 Catalyst Queen's Road Belfast BT3 9DT	Contact: Seamus Watson E-Mail: Seamus.W@arc-regulatory.com Website: www.arc-regulatory.com
Unit 1 Concourse building 2 Catalyst Queen's Road Belfast BT3 9DT	Contact: Seamus Watson E-Mail: Seamus.W@arc-regulatory.com Website: www.arc-regulatory.com		

DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS and TISSUES	<u>Immunohistochemistry examination activities</u>	In-house documented procedures using manual methods or analysers in conjunction with manufacturers' instructions as specified
FFPE glass slides	Immunohistochemistry to detect CLDN 18	Roche Ventana Benchmark UltraPlus Roche CLDN 18.2 kit Lab Working instruction 01834 IHC-CLDN18 protocol
FFPE glass slides	Digital imaging of stained-glass slides	Roche Ventana DP 200 Digital Scanner, Ventana Image Viewer and Roche Navify Digital Pathology Software Lab Work Instruction-00930
Digital images	Morphological assessment and interpretation/diagnosis	Digital Pathology SOP-00935 Reporting of results
Plasma extracted from Whole Blood in K2EDTA tubes	Molecular examination activities Molecular detection of EGFR mutation	Roche Cobas z480 Roche Cobas EGFR mutation kit Lab Working Instruction 01783 Cobas EGFR Mutation test
Serum extracted from whole blood	Serological examination activities Detection of AAV9 antibodies	TECAN infinite Nano and Plate reader and the ACRO Biosystems AAV9 ELISA kit Lab Work Instruction 01778 AAV9 ELISA
END		