


Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 <p>UKAS MEDICAL 8067</p> <p>Accredited to ISO 15189:2022</p>	<p>North Bristol NHS Trust</p> <p>Issue No: 013 Issue date: 17 February 2026</p>	
	<p>Department of Immunology and Immunogenetics Southmead Hospital Southmead Road Westbury on Trym Bristol BS10 5NB</p>	<p>Contact: Charlotte Besant Tel: +44 (0)117 4148372 E-Mail: charlotte.besant@nbt.nhs.uk Website: www.nbt.nhs.uk</p>
<p>Testing performed at the above address only</p>		

DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS</p> <p>Serum unless otherwise specified</p>	<p><u>Immunology examination activities for the purpose of clinical diagnosis</u></p> <p>Detection and quantification of specific IgE to the following allergens:</p> <p>Almond Apple rArah2 Aspergillus fumigatus IgE Banana Birch Blue Mussel Brazil Nut Cashew Nut Cat Epithelium Cheese, cheddar Cheese, mould Chilli Pepper Chocolate (Cacao) Coconut Codfish Dog Dander Egg White Gluten Grapefruit Guinea Pig Epithelium Hamster Epithelium Hazel Nut</p>	<p>Manufacturer's instructions for Thermofisher Immunocap 250 analyser (fluorescence solid phase immunoassay) and in-house procedure:</p> <p>IMMALSOP003</p>



8067

Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS (cont'd)</p> <p>Serum unless otherwise specified (cont'd)</p>	<p><u>Immunology examination activities for the purpose of clinical diagnosis</u> (cont'd)</p> <p>Detection and quantification of specific IgE to the following allergens: (cont'd)</p> <p>Honey Bee Horse Dander House Dust Mite Kiwi Fruit Latex Lemon Milk Mouse Epithelium Oak Orange Ovomucoid Peanut Pecan Nut Penicilloyl G Penicilloyl V Pistachio Potato Rabbit Epithelium Rat Epithelium Salmon Sesame Seed Shrimp (Prawn) Soya Bean Strawberry Sycamore Timothy Grass Tomato Tuna Walnut Wasp Wheat White Bean Fungal Antigen Aspergillus IgG</p>	<p>Manufacturer's instructions for Thermofisher Immunocap 250 analyser fluorescence solid phase immunoassay and in house procedure:</p> <p>IMMALSOP003</p>



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS (cont'd)</p> <p>Serum unless otherwise specified (cont'd)</p>	<p><u>Immunology examination activities for the purpose of clinical diagnosis</u> (cont'd)</p> <p>Detection and quantification of specific IgE to allergens - mixed antigen panels</p> <p>Mixed Nuts panel (Peanut, Hazel nut, Brazil nut, Almond, Coconut)</p> <p>Mixed Feathers panel (Goose, Chicken, Duck, Turkey)</p> <p>Mixed Cage Birds panel (Budgerigar, Parrot, Canary, Parakeet, Finch)</p> <p>Mixed Seafood panel (Cod, Shrimp, Blue mussel, Tuna, Salmon)</p> <p>Mixed Cereals panel (Wheat, Oat, Maize, Sesame seed, Buckwheat)</p> <p>Mixed Foods panel (Egg white, milk, cod, wheat, peanut, soya)</p> <p>Mixed Meats panel (Pork, Beef, Chicken and Turkey)</p> <p>Mixed Moulds(Penicillium, Cladosporium, Aspergillus, and Alternaria)</p> <p>Mixed Weeds (Mugwort, Plantain, Goosefoot, Golden Rod, and Nettle)</p> <p>Detection and quantification of Total IgE</p> <p>Mast Cell Tryptase</p>	<p>Manufacturer's instructions for Thermofisher Immucap 250 analyser fluorescence solid phase immunoassay and in-house procedure:</p> <p>IMMALSOP003</p>



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd) Serum unless otherwise specified (cont'd)	<u>Immunology examination activities for the purpose of clinical diagnosis</u> (cont'd)	Manufacturer's instructions for Thermofisher Immunocap 250 analyser (ELiA) and in-house procedure:
	Detection of antibodies to dsDNA (quantitative)	IMMAISOP028
	Detection of Anti-Transglutaminase (tTG) antibodies	IMMAISOP028
	Detection of Thyroids Stimulating Hormone Receptor Antibodies	IMMAISOP028
		Manufacturer's instructions for Quantalyser stainer and fluorescence microscopy or manual staining assay with fluorescence microscopy and in-house procedure:
	Detection of Anti-neutrophil Cytoplasmic Antibodies (ANCA)	IMMAISOP030 IMMAISOP006 (manual ANCA)
	ANCA titre	IMMAISOP030 IMMAISOP006 (manual ANCA)
	Auto Immune Profile (AIP)	IMMAISOP030 IMMAISOP005 (manual AIP)
	Detection of Anti-Nuclear Antibodies (Hep2 cells)	IMMAISOP030
	Detection of organ specific autoantibodies Adrenal Islet cell Paraneoplastic Skin IgA/IgG Endomysial Ds DNA using Crithidia	IMMAISOP033 IMMAISOP011 (manual)



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS (cont'd)</p> <p>Serum unless otherwise specified (cont'd)</p>	<p><u>Immunology examination activities for the purpose of clinical diagnosis (cont'd)</u></p> <p>Detection of Anti Cardiolipin IgG antibodies</p> <p>Detection of Anti-Cyclic Citrullinated Peptide (CCP) antibodies</p> <p>Detection of Anti PLA2-R antibodies</p>	<p>Manufacturer's instructions for Dynex DS2 (ELISA) and in-house procedure:</p> <p>IMMAISOP021</p> <p>IMMAISOP021</p> <p>IMMAISOP027</p>



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS (cont'd)</p> <p>Serum unless otherwise specified (cont'd)</p>	<p><u>Immunology examination activities for the purpose of clinical diagnosis (cont'd)</u></p> <p>Detection of liver antibodies M2 LKM-1 LC-1 and SLA/LP</p> <p>Detection of Paraneoplastic autoantibodies</p> <p>Detection of systemic sclerosis antibodies (Scl-70, CENP A, CENP B, RP11, RP155, fibrillarin, NOR90, Th/To, PM-Sc100, PM-Sc175, Ku,)</p> <p>Detection of Myositis antibodies (Mi-2 alpha, Mi-2 beta, TIF1g, MDA5, NXP2, SAE1, Ku, PM-Sc100, PM-Sc175, Jo-1, SRP, PL-7, PL-12, EJ, OJ, Ro-52)</p> <p>Detection of Extractable Nuclear Antigen specificities (ENA) :</p> <p>Anti Ro antibody Anti La antibody Anti Sm antibody Anti RNP antibody Anti Scl-70 antibody Anti Jo-1 antibody</p> <p>Determination of serum proteins:</p> <p>IgG Subclasses C3/C4 C1Inhibitor</p> <p>Free Light chains</p>	<p>Manufacturers instructions for Euroimmun Euroblot analyser (line blot) and manual method and in-house procedure:</p> <p>IMMAISOP035 (automated) IMMAISOP020 (manual)</p> <p>IIMMAISOP035 (automated) IMMAISOP018 (manual)</p> <p>IIMMAISOP035 (automated) IMMAISOP032 (manual)</p> <p>IIMMAISOP035 (automated) IMMAISOP031(manual)</p> <p>IMMAISOP035 (automated) IMMAISOP036 (manual)</p> <p>Manufacturer's instructions for Binding Site Optilite (turbidimetric assay) and in house procedure:</p> <p>IMMICSOP030</p> <p>IMMICSOP030</p>



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	<u>Immunology examination activities for the purpose of clinical diagnosis</u> (cont'd)	Manufacturers instructions for the Beckman Coulter AU 5800 in Biochemistry (UKAS 8071) and in house procedure:
Serum unless otherwise specified (cont'd)		
Serum and Plasma	Immunoglobulins: IgG, IgA and IgM	BS/CB/AUTO/BECK/12
Serum and Plasma	Beta-2 microglobulin	BS/CB/AUTO/BECK/12
Serum	Alpha-1 Antitrypsin	BS/CB/AUTO/BECK/12
		Manufacturers instructions for Fujirebio Lumipulse (chemiluminescence) analyser and in-house procedure
	Detection of Mesothelin antibodies	IMMAISOP034
		Manufacturer's instructions for DYNEX DS2 (ELISA) and in-house procedure:
Faeces	Detection of Faecal Elastase	IMMAISOP021 IMMICSOP019
Faeces	Detection of Faecal Calprotectin	IMMAISOP021 IMMICSOP018 using the Calex extraction device
		Manufacturers instructions for HM JACKarc Integrated Sphere Latex Turbidimetry and in-house procedure
Faeces	Quantitation of faecal haemoglobin	IMMICSOP033



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd) Serum unless otherwise specified (cont'd)	<u>Immunology examination activities for the purpose of clinical diagnosis</u> (cont'd)	Manufacturer's instructions for SEBIA CAPILLARYS 3 TERA(capillary electrophoresis) and in-house procedure:
	Detection of serum monoclonal proteins	IMMICSOP014
	Identification of serum monoclonal proteins	Manufacturer's instructions for SEBIA HYDRASYS (gel electrophoresis and immunofixation) and in-house procedure: Serum Immunofixation Electrophoresis IMMICSOP006
Urine	Detection of urine monoclonal proteins	Urine Immunofixation Electrophoresis (IFE) IMMICSOP021/IMMICSOP012
Serum, Body fluids	Detection of Asialotransferrin	IMMICSOP023
	Alpha 1 antitrypsin phenotype	IMMICSOP027
	C3 nephritic factor	IMMICSOP004
Whole Blood (EDTA)	C3 degradation products	IMMICSOP004
Serum	Cryoglobulins	Serum Immunofixation Electrophoresis (IFE) IMMICSOP005
Serum	Alkaline Phosphatase Isoenzyme	Sebia Hydrasys 2 Electrophoresis IMMIC SOP040
Serum and CSF	Detection of Oligoclonal Bands	Sebia Hydrasys 2 Electrophoresis IMMIC SOP039



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	<u>Immunophenotyping examination activities for the purpose of clinical diagnosis</u>	Manufacturer's instructions for Becton Dickinson Facs Lyric (flow cytometry) and in-house procedure:
Whole Blood (EDTA), Bone Marrow, Lymph Node, Body Fluids	Leukaemia/lymphoma type	IMMIPSOP007 GP/PPP/0017
	CD1a	IMMIPSOP008
	CD2	IMMIPSOP012
	CD3	IMMIPSOP029
	CD44	IMMIPSOP043
	CD5	IMMIPSOP044
	CD7	
	CD8	
	CD10	
	CD11c	
	CD11b	
	CD13	
	CD14	
	CD15	
	CD19	
	CD20	
	CD22	
	CD23	
	CD25	
	CD26	
	CD30	
	CD33	
	CD34	
	CD36	
	CD38	
	CD41a	
	CD42b	
	CD43	
	CD45	
	CD45RA	
	CD56	
	CD57	
	CD64	
	CD71	
	CD79a	
	CD79b	
	CD103	
	CD117	
	CD123	
	CD200	
	cMPO	



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	<u>Immunophenotyping examination activities for the purpose of clinical diagnosis</u> (cont'd)	Manufacturer's instructions for Becton Dickinson Facs Lyric(flow cytometry) and in-house procedure:
Whole Blood (EDTA), Bone Marrow, Lymph Node, Body Fluids (cont'd)	Leukaemia/lymphoma type (cont'd)	IMMIPSOP007 GP/PPP/0017
	cTdt cIgM cCD3 DR GlycoA Kappa Lamda	IMMIPSOP008 IMMIPSOP012 IMMIPSOP029 IMMIPSOP043 IMMIPSOP044 IMMIPSOP046
Whole Blood (EDTA), Bone Marrow, Lymph Node, Body Fluids	Detection of malignant plasma cells CD19 CD27 CD38 CD45 CD56 CD81 CD117 CD138	IMMIPSOP009 IMMIPSOP007 GP/PPP/0017 IMMIPSOP008 IMMIPSOP012 IMMIPSOP029 IMMIPSOP043 IMMIPSOP044 IMMIPSOP046



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS (cont'd)</p> <p>Whole Blood (EDTA), Bone Marrow, Lymph Node, Body Fluids</p>	<p><u>Immunophenotyping examination activities for the purpose of clinical diagnosis</u> (cont'd)</p> <p>Detection of Minimal Residual Disease (AML, ALL, CLL)</p> <p>CD1a CD2 CD3 CD4CD5 CD7 CD8 CD10 CD13 CD15 CD16 CD19 CD20 CD22 CD33 CD34 CD38 CD45 CD48 CD56 CD58 CD66c CD73 CD79b CD81 CD99 CD117 CD123 CD304 NG2 Dr cCD3</p>	<p>Manufacturer's instructions for Becton Dickinson Facs Lyric(flow cytometry) and in-house procedure:</p> <p>IMMIPSOP021 IMMIPSOP007 GP/PPP/0017 IMMIPSOP008 IMMIPSOP012 IMMIPSOP043 IMMIPSOP044 IMMIPSOP046</p>



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	<u>Immunophenotyping examination activities for the purpose of clinical diagnosis (cont'd)</u>	Manufacturer's instructions for Becton Dickinson Facs Lyric(flow cytometry) and in-house procedure:
Whole Blood (EDTA)	Immunodeficiency screen CD3 CD4 CD8 CD19 CD45 CD56/16	IMMIPSOP045 IMMIPSOP036
Whole Blood (EDTA)	Rituximab monitoring CD19	IMMIPSOP036 IMMIPSOP045
Whole Blood (EDTA)	Quantification of memory B cells IgD IgM CD19 CD27	IMMIPSOP027 IMMIPSOP007 IMMIPSOP012 IMMIPSOP043 IMMIPSOP044
Whole Blood (EDTA)	Quantification of dual negative CD4/CD8 cells CD3 CD4 CD8 TCR (alpha/beta)	IMMIPSOP028 IMMIPSOP007 IMMIPSOP012 IMMIPSOP043 IMMIPSOP044
Whole Blood (EDTA)	Quantification of T cell subsets CD3 CD4 CD8 CD45	IMMIPSOP015 IMMIPSOP045
Whole Blood (EDTA)	Identification of PNH clones FLAER CD14 CD15 CD45 CD55 CD59 GlycoA	IMMIPSOP010 IMMIPSOP007 GP/PPP/0017 IMMIPSOP012 IMMIPSOP043 IMMIPSOP044



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	<u>Immunophenotyping examination activities for the purpose of clinical diagnosis (cont'd)</u>	Manufacturer's instructions for Becton Dickinson Facs Lyric(flow cytometry)and in-house procedure:
Whole Blood (Heparin)	Functional Assays Neutrophil function- Adhesion molecules	IMMIPSOP012 IMMIPSOP043 IMMIPSOP044 IMMIPSOP048
Whole Blood (Heparin)	Lymphocyte function tests	IMMIPSOP011
Whole Blood (Heparin)	CD40 Ligand	IMMIPSOP011
Whole Blood (Test specific samples vials)	Tb interferon gamma release assay	Manufacturer's instructions for DYNEX DS2 (ELISA) and QUANTIFERON assay and in-house procedure: IMMIPSOP041



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
<p>HUMAN TISSUE AND FLUIDS (cont'd)</p> <p>Serum unless otherwise specified (cont'd)</p> <p>Whole Blood (EDTA or Citrate), Lymph Node, Spleen</p>	<p><u>Histocompatibility and Immunogenetics examination activities for the purpose of clinical diagnosis</u></p> <p>DNA Extraction</p>	<p>Promega Maxwell RSC automated extraction platform in-house procedure:</p> <p>TTTYP SOP046</p>
<p>DNA. extracted from sample types listed above</p>	<p>PCR HLA typing HLA class I/II HLA-A, B, C, DRB1, DRB3/4/5, DQB1, DQA1, DPB1 and DPA1</p> <p>HLA B57/B*57:01</p> <p>HLA B27 Ankylosing Spondylitis</p> <p>HLA B51 (5) Behcets Syndrome</p> <p>HLA DR15(2), DQ6(DQB1*06:02) Narcolepsy</p> <p>HLA A29 Birdshot Retinopathy</p> <p>HLA-DQ2 and HLA DQ8 Coeliac disease</p>	<p>Manufacturer's instructions for PCR and Luminex FlexMap 3D PCR-SSO assay using Lifecodes kit and MatchIt software</p> <p>TTTYP SOP 048 TTTYP SOP049</p>
<p>DNA. extracted from sample types listed above</p>	<p>PCR HLA typing HLA class I/II HLA-A, B, C, DRB1, DRB3/4/5, DQB1, DQA1, DPB1 and DPA1</p>	<p>Manufacturer's instructions for Roche Lightcycler 480 II RT-PCR and CareDX QType assay and in-house procedure</p> <p>TTTYP SOP050</p>



8067
Accredited to
ISO 15189:2022

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

North Bristol NHS Trust
Issue No: 013 Issue date: 17 February 2026

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	<u>Histocompatibility and Immunogenetics examination activities for the purpose of clinical diagnosis (cont'd)</u>	Manufacturer's instructions for Luminex FlexMap 3D (multiplex bead assay with fluorescence detection) using One Lambda kits HLA Fusion software
Serum unless otherwise specified (cont'd)		
Serum	HLA antibody screen	TTSCRSOP005 TTSCRSOP006
Serum	HLA antibody identification (Class I and Class II single antigen beads)	TTSCRSOP005
Serum	HLA antibody identification (Class I and Class II single antigen beads)	TTSCRSOP006
Serum	HLA incompatibility transplantation (Class I and Class II single antigen beads)	TTSCRSOP016
Whole Blood (EDTA or Citrate) Lymph node, spleen	Flow Cytometry Crossmatch	Becton Dickinson Facs Lyric (Flow cytometry) TTRENSOP008 TTRENSOP017
Whole Blood (EDTA)	Monitoring of Desensitization to ABO incompatible Donor	Manual In house procedure: (Manual titration) TTRENSOP016

END