Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK



9595

Accredited to ISO 15189:2022

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Chester & Wirral Microbiology

Service

11 Bassendale Road

Wirral CH62 3QL Contact: Mr Alex Warrington Tel: +44 (0)1244 362499

E-Mail: alex.warrington@nhs.net

Testing performed at the above address only

DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUES AND FLUIDS	Microbiological examination activities for the purpose of clinical diagnosis	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Routine samples: Aspirates Fluids Pus Samples from sterile sites Urine Bone Cellular and fluid material contained on swabs, catheter tips and intra-uterine contraceptive devices Genital samples Semen Tissues	Isolation of microorganisms of clinical significance	Manual inoculation and media culture using: MBSOP 33 MBSOP 21 MBSOP 18 MBSOP 29
Cellular and fluid material contained on swabs	Screening for Methicillin- Resistant <i>Staphylococcus</i> <i>aureus</i> (MRSA)	Manual inoculation and media culture using: MBSOP 44
Rectal swab, Faeces and Urine	Screening for Resistant organisms VRE and CPE	Manual inoculation and media culture using: MBSOP 46
Blood Joint fluids	Isolation of microorganisms of clinical significance	Liquid media culture using: MBSOP 52 MESOP 12 and Biomérieux BacTalert



Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	Microbiological examination activities for the purpose of clinical diagnosis (cont'd)	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Cerebrospinal fluid (CSF)	Detection, identification and quantification (inc. differential)	Microscopy & staining (Gram and differential) using:
	of: White blood cells Red blood cells Bacterial cells	MBSOP 11
CSF	Isolation of microorganisms of clinical significance	Manual inoculation and media culture using:
		MBSOP 11
Bronchoalveolar lavage Bronchial aspirate, brushings & washings Cellular and fluid material contained on swabs Endotrachael tube specimens Nasopharyngeal aspirate Pleural fluid Sputum	Isolation of microorganisms of clinical significance	Manual inoculation and media culture using: MBSOP 42
	Detection, identification and quantification of:	Staining (inc. Gram) and microscopic examination using:
Aspirates Fluids CAPD Corneal scrapes Eye (specimens from) Tissues BAL	White blood cells Red blood cells Bacterial cells Yeast cells and fungal hyphae	MBSOP 21 MBSOP52 MBSOP42



Accredited to ISO 15189:2022

Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	Microbiological examination activities for the purpose of clinical diagnosis (cont'd)	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Urine	Detection, identification and quantification of: White blood cells Red blood cells Squamous epithelial cells Bacterial cells Yeasts	Automated urinalysis, and manual examination using: MBSOP 33 and SediMAX contrast (x2) MESOP 50 and/or Manual microscopic examination using: KOVA slides
Urine	Isolation of microorganisms of clinical significance	Manual inoculation and media culture using: MBSOP 33 BG9
Faeces	Isolation of bacterial pathogens: Campylobacter spp. Escherichia coli O157 Salmonella sp. Shigella spp. Vibrio spp. Yersinia spp. Pleisiomonas spp.	Manual inoculation and media culture using: MBSOP59
Hair Nail Skin	Isolation and identification of clinically significant fungi	Manual inoculation, media culture, and macroscopic & microscopic examination using: MBSOP 38
Cultures isolated in-house from all samples listed above	Identification of microorganisms of clinical significance	Manual and/or automated biochemical methods using: Vitek 2- MESOP 24 and Matrix-assisted laser desorption/ionisation using: MESOP 32 Vitek MS Maldi-ToF



ISO 15189:2022

Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	Microbiological examination activities for the purpose of clinical diagnosis (cont'd)	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Bacterial cultures isolated inhouse from all samples listed above	Antimicrobial susceptibility testing of clinically significant bacteria	Manual and/or automated antimicrobial susceptibility testing using the latest Global EUCAST methodology Vitek 2 using in-house procedure MESOP 24 EUCAST Disc diffusion & gradient minimum inhibitory concentration determination using inhouse procedure MBSOP 06 and DOC 2689
Female genital tract swabs	Microscopic diagnosis of bacterial vaginosis	Gram stained microscopic examination using: MBSOP 29 Automated staining using Prolab Polystainer as per MESOP39
Faeces	Detection of ova, cysts and parasites (OCP)	Unstained and iodine stained direct microscopic examination and/or Formol-ether faecal concentration and microscopic examination using: MBSOP 59
Faeces	Detection of: Giardia sp Cryptosporidium sp Campylobacter sp Salmonella enterica Shigella sp. E.coli (STEC) C.difficile Toxin B detection	Enteric PCR using: MBSOP78, MBSOP59 and Serosep EntericBio Workstation Roche Lightcycler 480 II
Urine	Detection of Schistosoma haematobium ova	Microscopic examination using in-house procedure: MBSOP 33



Accredited to ISO 15189:2022

Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUE AND FLUIDS (cont'd)	Microbiological examination activities for the purpose of clinical diagnosis (cont'd)	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Nose and throat Swabs Faeces Culture isolate suspension	Detection or exclusion of: Influenza A Influenza B Respiratory syncytial virus MRSA CPE VRE Norovirus Difficile Toxin Genes tcdB, cdt and tcdC	Automated Molecular Testing (PCR) using: MBSOP 73 MBSOP 49 MBSOP 71 MBSOP 72 MVSOP 14 MBSOP 70 and: Cepheid GeneXpert
Nose/throat swabs in VTM	Detection of SARS-CoV-2 virus RNA (COVID-19) specific gene sequences: E, N2	Cepheid GeneXpert MBSOP79 Cepheid PCR Xpert® Xpress SARS-CoV-2 (Ref XPRSARS-COV2-10)
Nose/throat swabs Nasal swabs	Detection of: SARS-CoV-2(RdRP, N2 & E gene targets) Influenza A Influenza B RSV	Cepheid GeneXpert (PCR) Xpert Xpress combined Multiplex assay MBSOP84
Urine	Detection of Legionella pneumophilia serotype 1 antigen and Streptococcus pneumoniae antigen	Immuno-chromatographic assay using: MVSOP 12 and Launch TRU Urine Antigen kit
Urine, endocervical/vaginal swabs, rectum and throat	Detection or exclusion of: Chlamydia trachomatis and/or Neisseria gonorrhoeae	Multiplex PCR using MBSOP7 and Roche Cobas 4800
Faeces	Rotavirus and adenovirus	Immuno-chromatographic assay Coris Combi-strip MVSOP 13
Faeces	Detection of: GDH Toxins A & B	Techlab C.difficile Quik Chek Complete kit Rapid membrane enzyme immunoassay MBSOP70



ISO 15189:2022

Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUES AND FLUIDS (cont'd)	Serological examination activities for the purpose of clinical diagnosis	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Faeces	Helicobacter Antigen	EIA – Dynex DS2 Premier Platinum HpSA Plus MBSOP 09 MESOP 29 ** Documented in house methods to meet the requirements of the Infectious Diseases in pregnancy Screening (IDPS) as defined in the July 2018 'NHS Infectious Diseases in Pregnancy Programe: laboratory quality assurance evidence requirements'
Plasma Serum	Microbial serology for: Hepatitis B surface antibody CMV IgM and IgG & avidity Hepatitis A IgG and IgM **Hepatitis B surface antigen Hepatitis C antibody **HIV-1 & HIV-2 Antibodies and HIV-1 p24 Antigen HTLV ab Rubella IgG **Syphilis Total Antibody (IgG & IgM) Toxoplasma IgG Hepatitis B Core total HBV markers	Enzyme Immuno-Assay (EIA) Roche Cobas 6000 using: MESOP 30 MVSOP 29 MVSOP 34 MVSOP 04 MVSOP 02 MVSOP 28 MVSOP 05 MVSOP 27 MVSOP 07 MVSOP 15 MVSOP 09 MVSOP29



ISO 15189:2022

Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUES AND FLUIDS (cont'd)	Serological examination activities for the purpose of clinical diagnosis (cont'd)	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Plasma Serum	Microbial serology for: Chlamydia IgG EBV IgM and IgG Measles IgG Varicella Zoster IgG Mycoplasma IgM Parvovirus IgM and IgG EBV (Epstein Barr Virus): EBV VCA IgG EBV VCA IgM EBNA IgG Mycoplasma IgM	EIA - Dynex BEST 2000 using: MESOP 09 MVSOP 32 MVSOP 25 MVSOP 30 MVSOP 06 MVSOP 10 DOC 1613 MVSOP25
Plasma Serum	Microbial serology for:	EIA - VIDAS using: DOC1619
	Hepatitis B Markers (HBe, Ab and HBsAg) Varicella IgG Hepatitis B Core IgM Lymes (Borrelia) IgM and IgG Measles IgG Hepatitis C Ab HIV 1&2	MVSOP 02 MVSOP 06 MVSOP 02 MVSOP 24 MVSOP 30 MVSOP 28 MVSOP 05
Plasma Serum	HIV antigen-antibody	Immuno-chromatographic strip method Alere HIV-combo kit MVSOP 05
Serum/plasma	Syphilis antibody confirmation	Abbott Bioline Manual lateral flow MVSOP 15



Schedule of Accreditation issued by

United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Wirral University Teaching Hospital NHS Foundation Trust

Issue No: 009 Issue date: 11 November 2025

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN TISSUES AND FLUIDS (cont'd)	Serological examination activities for the purpose of clinical diagnosis (cont'd)	In-house documented procedures following current UK Standards for Microbiology Investigations (UK SMIs)
Nose & Throat swabs Nasal swabs	Adenovirus Coronavirus 229E Coronavirus HKU1 Cornavirus NL63 Coronavirus OC43 Middle East Respiratory Syndrome Coronavirus (MERS-CoV) Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) Human Metapneumovirus Human Rhinovirus/Enterovirus Influenza A, including subtypes H1, H3 & H1-2009 Influenza B Parainfluenza Virus 1 Parainfluenza Virus 2 Parainfluenza Virus 3 Parainfluenza Virus 4 Respiratory Syncytial Virus (RSV) Bordetella parapertussis Bordetella pertussis Chlamydia pneumoniae Mycoplasma pneumonia	Biofire filmArray Multiplex PCR MBSOP88
END		