


# Schedule of Accreditation

issued by

## United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 <b>9766</b> <b>Accredited to ISO 15189:2012</b>	<b>Western Health and Social Care Trust</b>	
	<b>Issue No:</b> 004 <b>Issue date:</b> 18 February 2025	
	<b>Department of Haematology and Blood Transfusion Laboratory South West Acute Hospital  Irvinestown Road Enniskillen BT74 4RT</b>	<b>Contact: Adrian Crawford</b> <b>Tel: +44 (0) 2871 345 171 ext 213789</b> <b>Fax: +44 (0) 2871 313036</b> <b>E-Mail: <a href="mailto:Adrian.crawford@westerntrust.hscni.net">Adrian.crawford@westerntrust.hscni.net</a></b>
Testing performed at the above address only		

### Site activities performed away from the location listed above:

Location details		Activity
<b>Address</b>  Blood Transfusion Laboratory, South West Acute Hospital, Enniskillen	<b>Local Contact</b>  Adrian Crawford	Management of blood fridges  3 Blood Fridges, 2 Freezers
Tyrone County Hospital Cardiac Assessment Unit Tyrone County Hospital Omagh	Adrian Crawford	Management of blood fridges  Cardiac Assessment Unit (1 fridge)



9766

Accredited to  
ISO 15189:2012

## Schedule of Accreditation

issued by

### United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Western Health and Social Care Trust

Issue No: 004 Issue date: 18 February 2025

Testing performed at main address only

### DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS	<u>Haematology examination activities for the purpose of clinical diagnosis</u>	Manufacturer's instructions and documented In-house methods using:
Blood	Full Blood Count <ul style="list-style-type: none"> <li>- Haemoglobin</li> <li>- White blood cell count</li> <li>- Red blood cell count</li> <li>- Platelet count</li> <li>- White blood cell differential count</li> <li>- Reticulocyte count</li> <li>- Nucleated Red Blood Cell Count (NRBC)</li> <li>- Mean corpuscular volume</li> <li>- Mean corpuscular haemoglobin</li> <li>- Mean corpuscular haemoglobin concentration</li> <li>- Haematocrit</li> <li>- Red cell distribution width</li> </ul>	Sysmex XN-20 analyser  Key SOPs Haem/3412Haem/18 -> 42
Blood films	Reporting of normal and abnormal morphologies	Stainers: Siemens Hematek stainer  Haematology Microscopes: Nikon E600 & Olympus BX51.  Key SOPs Haem/33-36 Haem 42 Haem/290-295
Blood	Erythrocyte Sedimentation Rate (ESR)	StaRRsed Interrliner ESR analyser.  Principle: Westergren method  Key SOPs Haem/37-40, Haem/296
Blood	Screening for the presence of malaria parasites	Carestart Malaria RDT Kit.  Thin films using manual Giemsa stain  Key SOPs Haem/299 Haem/300



9766  
Accredited to  
ISO 15189:2012

**Schedule of Accreditation**  
issued by  
**United Kingdom Accreditation Service**  
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

**Western Health and Social Care Trust**

**Issue No: 004    Issue date: 18 February 2025**

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	<u>Haematology examination activities for the purpose of clinical diagnosis</u> (cont'd)	Manufacturer's instructions and documented In-house methods using:
Blood	Detection of infectious mononucleosis heterophile antibodies	Clearview Infectious Mononucleosis (IM) Kit.  Qualitative immunoassay kit.  Key SOP Haem/298
Blood	Detection of haemoglobin S (Hb-S),	Sickle-chex kit.  Principle: relative insolubility of Haemoglobin S when combined with buffer and a reducing agent  Key SOP Haem/ 303
HUMAN BODY FLUIDS	<u>Blood Transfusion examination activities for the purpose of clinical diagnosis</u>	Manufacturer's instructions and documented In-house methods using:
Blood	Blood Grouping  A, B, AB, O	Automated methods using Bio-Rad IH-1000 Analysers Manual Bio-Rad Gelcard method Manual Tube Method  Key SOPs BloodBank 135 Bloodbank 104 Bloodbank/37/38/39



9766

Accredited to  
ISO 15189:2012**Schedule of Accreditation**

issued by

**United Kingdom Accreditation Service**

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

**Western Health and Social Care Trust****Issue No:** 004 **Issue date:** 18 February 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	<u>Blood Transfusion examination activities for the purpose of clinical diagnosis (cont'd)</u>	Manufacturer's instructions and documented In-house methods using:
Blood	Screening for clinically significant antibodies of the following systems: <ul style="list-style-type: none"><li>• Rhesus</li><li>• Kell</li><li>• Duffy</li><li>• Kidd</li><li>• Lewis</li><li>• MNS</li><li>• P</li><li>• Lutheran systems</li></ul>	Automated methods using Bio-Rad IH-1000 Analysers Manual Bio-Rad Gelcard method  Key SOPs BloodBank 135 Bloodbank 77
Blood	Antigen phenotyping and genotyping of clinically significant antigens of the following systems: <ul style="list-style-type: none"><li>• Rhesus</li><li>• Kell</li><li>• Duffy</li><li>• Kidd</li><li>• Lewis</li><li>• MNS</li><li>• P</li><li>• Lutheran systems</li></ul>	Key SOPs Bloodbank 135 Bloodbank 83
Blood	Group & Coombs for blood grouping of babies up to 4 months: A, B, AB, O	Key SOPs BloodBank 135 Bloodbank 104
Blood	Cross matching	Methods as above plus Manual Tube method.  Key SOPs BloodBank 135 Bloodbank 78



9766  
Accredited to  
ISO 15189:2012

**Schedule of Accreditation**  
issued by  
**United Kingdom Accreditation Service**  
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

**Western Health and Social Care Trust**

**Issue No: 004    Issue date: 18 February 2025**

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	<u>Blood Transfusion examination activities for the purpose of clinical diagnosis (cont'd)</u>	Manufacturer's instructions and documented In-house methods using:
Blood	Direct Antiglobulin test	Automated methods using Bio-Rad IH-1000 Analysers Manual Bio-Rad Gelcard method Key SOPs  BloodBank 135 Bloodbank 43
Blood	Indirect Antiglobulin test	Key SOPs  BloodBank 135 Bloodbank 43
Blood	Identification of clinically significant antibodies of the following systems: <ul style="list-style-type: none"> <li>• Rhesus</li> <li>• Kell</li> <li>• Duffy</li> <li>• Kidd</li> <li>• Lewis</li> <li>• MNS</li> <li>• P</li> <li>• Lutheran systems</li> </ul>	Key SOPs  BloodBank 79 BloodBank 135 Bloodbank 79
Blood	Kleihauer Test for estimation of foeto-maternal haemorrhage	Manual staining method using Inverclyde Foetal Red Cell Detection Kit.  Key SOP BloodBank1-135



9766

Accredited to  
ISO 15189:2012

## Schedule of Accreditation

issued by

### United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Western Health and Social Care Trust

Issue No: 004 Issue date: 18 February 2025

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS (cont'd)	<u>Coagulation examination activities for the purpose of clinical diagnosis</u>	
Blood	Blood Coagulation: <ul style="list-style-type: none"><li>• PT</li><li>• APTT</li><li>• Fibrinogen</li><li>• D-Dimer</li><li>• International Normalised Ratio (INR)</li><li>• </li></ul>	Sysmex ACL TOP 550 analyser  Key SOPs  Coag/1008 Coag/327 - 334
Urine	<u>Miscellaneous examination activities for the purpose of clinical diagnosis</u>  Qualitative detection of human chorionic gonadotropin	Manufacturer's instructions and documented In-house methods using:  Alere Pregnancy Test Kit.  Key SOP Haem/297
END		