## **Schedule of Accreditation**

issued by

## **United Kingdom Accreditation Service**

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK



Randox Laboratories Limited, Trading as Randox Testing Services

Issue No: 039 Issue date: 11 July 2023

4348

Accredited to ISO/IEC 17025:2017

**BT29 4QY** 

 34 Diamond Road
 Contact: Mr Niall Fox

 Crumlin
 Tel: +44 (0)2894 422413

 Co Antrim
 Fax: +44 (0)1617 412777

E-Mail: Niall.Fox@randox.com

Website: www.randoxtestingservices.com

Testing performed at the above address only

#### **DETAIL OF ACCREDITATION**

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS	Workplace Drug Analysis and Medical Legal Analysis	
Urine (Preserved and unpreserved)	Presumptive Screening for the presence of drugs of abuse (DOAI and DOA II) (cut-off limit): Amphetamines (500 ng/mL) Methamphetamine (500 ng/mL) Benzodiazepines (200 ng/mL) Barbiturates (200 ng/mL) Opiates (300 ng/mL) Cocaine metabolites (150 ng/mL) Methadone (300 ng/mL) Cannabinoids (50 ng/mL) Phencyclidine (25 ng/mL) Ketamine (75ng/mI) Oxycodone Class 1 and Class 2 (100ng/mI) Generic opioids (100ng/mI) FentanyI (1ng/mI) Propoxyphene (300ng/mI) Methaqualone (300ng/mI)	Documented in house method SOP RTL 2213 RTL using Randox Evidence Investigator (method reference EI) biochip assay

Assessment Manager: FB Page 1 of 9



Accredited to ISO/IEC 17025:2017

### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, **Trading as Randox Testing Services**

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Urine (Preserved and unpreserved)	Presumptive Screening for the presence of drugs of abuse (cont'd) (DOAI and DOA II) (cut-off limit): Buprenorphine (2ng/ml) LSD (0.5ng/ml) MDMA Group (500ng/ml)	Documented in house method SOP RTL 2213 RTL using Randox Evidence Investigator (method reference EI) biochip assay
Urine (Preserved and unpreserved)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	
	Δ9-Tetrahydrocannabinoic acid (THC-COOH) (15ng/ml) [5-150ng/ml]	Documented in house method A4URN and WI0227-RTL using solid phase extraction LC-MS/MS
Urine (Preserved and Unreserved)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:  Amphetamines:-	Documented in house method WI1099-RTL using enzymatic hydrolysis (where applicable) and LC-MS/MS
	Amphetamine (200ng/ml) [30-800ng/ml] Methamphetamine (200ng/ml) [30-800ng/ml] MDMA (200ng/ml) [30-800ng/ml] MDA (200ng/ml) [30-800ng/ml]	

Assessment Manager: FB Page 2 of 9



## Accredited to ISO/IEC 17025:2017

### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, **Trading as Randox Testing Services**

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Urine (Preserved and Unreserved)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	Documented in house method WI1099-RTL using enzymatic hydrolysis (where applicable) and LC-MS/MS
	Cocainics:- Benzoylecgonine (100ng/ml) [30-800ng/ml] Cocaine (100ng/ml) [30-800ng/ml]	
	Opiates:- 6-monoacetylmorphine (10ng/ml) [3-80ng/ml] Morphine (300ng/ml) [30-800ng/ml] Codeine (300ng/ml) [30-800ng/ml] Dihydrocodeine (300ng/ml) [30-800ng/ml]	
	Opioids:- Methadone (250ng/ml) [30-800ng/ml] EDDP (75ng/ml) [30-800ng/ml] Oxycodone (300ng/ml) [30-800ng/ml] Oxymorphone (300ng/ml) [30-800ng/ml]	

Assessment Manager: FB Page 3 of 9



Accredited to ISO/IEC 17025:2017

## Schedule of Accreditation issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, Trading as Randox Testing Services

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Urine (Preserved and Unreserved)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	Documented in house method WI1099-RTL using enzymatic hydrolysis (where applicable) and LC-MS/MS
	Opioids (cont'd):- Hydrocodone (300ng/ml) [30-800ng/ml] Hydromorphone (300ng/ml) [30-800ng/ml] Fentanyl (1ng/ml) [0.3-8.0ng/ml] Tramadol (100ng/ml) [30-800ng/ml]  z-Drugs:- Zolpidem (10ng/ml) [3-80ng/ml]  Ketamine:- Ketamine (75ng/ml) [30-800ng/ml] Norketamine (75ng/ml) [30-800ng/ml]	

Assessment Manager: FB Page 4 of 9



## Accredited to ISO/IEC 17025:2017

### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, Trading as Randox Testing Services

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Urine (Preserved and Unreserved)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	Documented in house method WI1099-RTL using enzymatic hydrolysis (where applicable) and LC-MS/MS
	Benzodiazepines:- Alpha-hydroxyalprazolam (100ng/ml) [30-800ng/ml] 7-aminoclonazepam (100ng/ml) [30-800ng/ml] 7-aminoflunitrazepam (100ng/ml) [30-800ng/ml] 7-aminonitrazepam (100ng/ml) [30-800ng/ml] Bromazepam (100ng/ml) [30-800ng/ml] Desalkylflurazepam (100ng/ml) [30-800ng/ml] Diazepam (100ng/ml) [30-800ng/ml] Lorazepam (100ng/ml) [30-800ng/ml] Lormetazepam (100ng/ml) [30-800ng/ml] Nordiazepam (100ng/ml) [30-800ng/ml] Nordiazepam (100ng/ml) [30-800ng/ml] Oxazepam (100ng/ml) [30-800ng/ml] Oxazepam (100ng/ml) [30-800ng/ml] Phenazepam (100ng/ml) [30-800ng/ml] Phenazepam (100ng/ml) [30-800ng/ml]	

Assessment Manager: FB Page 5 of 9



Accredited to ISO/IEC 17025:2017

## Schedule of Accreditation issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, Trading as Randox Testing Services

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Urine (Preserved and Unreserved)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	Documented in house method WI1099-RTL using enzymatic hydrolysis (where applicable) and LC-MS/MS
	Benzodiazepines (cont'd):- Temazepam (100ng/ml) [30-800ng/ml]	
	Gabapentinoids:- Pregabalin (100ng/ml) [30-800ng/ml] Gabapentin (100ng/ml) [30-800ng/ml]	
	Ergoline:- LSD (1ng/ml) [0.3-8.0ng/ml]	
	Phencyclidine:- PCP (25ng/ml) [3-80ng/ml]	
	Quinazoline:- Methaqualone (100ng/ml) [30-800ng/ml]	

Assessment Manager: FB Page 6 of 9



Accredited to ISO/IEC 17025:2017

### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, **Trading as Randox Testing Services**

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS and TISSUES	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Human Hair (Head and Body)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	Documented in house method HDT and WI1790-RTL using LC- MS/MS
	Amphetamine (0.2ng/mg) [0.1-10ng/mg] Methamphetamine (0.2ng/mg) [0.1-10ng/mg] MDMA (0.2ng/mg) [0.1-10ng/mg] MDA (0.2ng/mg) [0.1-10ng/mg] Cocaine (0.5ng/mg) [0.25-25ng/mg] Norcocaine (0.05ng/mg) [0.025-2.5ng/mg] Cocaethylene (0.05ng/mg) [0.025-2.5ng/mg] m-Hydroxycocaine (0.05ng/mg) [0.025-2.5ng/mg] BZE (0.05ng/mg) [0.025-2.5ng/mg] AEME (0.05ng/mg) [0.025-2.5ng/mg] AEME (0.05ng/mg) [0.025-2.5ng/mg] Aeme (0.1-10ng/mg) [0.2ng/mg] 6-MAM (0.2ng/mg) [0.1-10ng/mg]	

Assessment Manager: FB Page 7 of 9



Accredited to ISO/IEC 17025:2017

### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, Trading as Randox Testing Services

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS and TISSUES (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Human Hair (Head and Body)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:	Documented in house method HDT and WI1790-RTL using LC- MS/MS
	Codeine (0.2ng/mg) [0.1-10ng/mg] Dihydrocodeine (0.2ng/mg) [0.1-10ng/mg] Alprazolam (0.05ng/mg) [0.025-2.5ng/mg] Tramadol (0.2ng/mg) [0.1-10ng/mg] Oxazepam (0.05ng/mg) [0.025-1.0ng/mg] Bromazepam (0.05ng/mg) [0.025-1.0ng/mg] Clonazepam (0.05ng/mg) [0.025-1.0ng/mg] Diazepam (0.05ng/mg) [0.025-1.0ng/mg] Nordiazepam (0.05ng/mg) [0.025-1.0ng/mg] Nordiazepam (0.05ng/mg) [0.025-1.0ng/mg] Nordiazepam (0.05ng/mg) [0.025-1.0ng/mg] Nordiazepam (0.05ng/mg) [0.025-1.0ng/mg] Phenazepam (0.05ng/mg) [0.025-1.0ng/mg] Phenazepam (0.05ng/mg) [0.025-1.0ng/mg] Phenazepam (0.05ng/mg) [0.025-1.0ng/mg]	

Assessment Manager: FB Page 8 of 9



Accredited to ISO/IEC 17025:2017

### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Randox Laboratories Limited, Trading as Randox Testing Services

Issue No: 039 Issue date: 11 July 2023

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
BODY FLUIDS and TISSUES (cont'd)	Workplace Drug Analysis and Medical Legal Analysis (cont'd)	
Human Hair (Head and Body)	Confirmation (above specified cut-off limit) and Quantitative analysis of the following drugs (cut-off limit) [concentration range]:  Temazepam (0.05ng/mg) [0.025-1.0ng/mg] Chlordiazepoxide (0.05ng/mg) [0.025-1.0ng/mg] Methadone (0.2ng/mg) [0.1-10ng/mg] EDDP (0.05ng/mg) [0.025-1.0ng/mg] Ketamine (0.5ng/mg) [0.25-25.0ng/mg] Norketamine (0.1ng/mg) [0.05-2.0ng/mg] THC (0.05ng/mg) [0.05-2.0ng/mg] THC (0.05ng/mg) [0.025-1.0ng/mg]	Documented in house method HDT and WI1790-RTL using LC-MS/MS
END		

Assessment Manager: FB Page 9 of 9