

Schedule of Accreditation

issued by

United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 <p>UKAS TESTING</p> <p>7774</p> <p>Accredited to ISO/IEC 17025:2017</p>	<p>Bourn Hall Limited, trading as Bourn Hall Clinic</p> <p>Issue No: 007 Issue date: 17 July 2020</p>	
	<p>Bourn Hall High Street Bourn CB23 2TN United Kingdom</p>	<p>Contact: Natalie Lowden Tel: 01954 719 111 Fax: 01954 718 826 E-Mail: info@bourn-hall.com Website: www.bourn-hall-clinic.co.uk</p>
<p>Testing performed at the above address only</p>		

DETAIL OF ACCREDITATION

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS	<p>Examination activities for the purposes of clinical diagnosis:</p> <p>Endocrinology</p> <p>Hormone analysis: quantification of: Anti-Mullerian hormone (AMH) beta-Human Chorionic Gonadotrophin (βHCG) Estradiol (E2) Follicle Stimulating Hormone (FSH) Leutenising Hormone (LH) Progesterone (P4) Testosterone (TEST) Thyroid stimulated hormone (TSH)</p> <p>Serology</p> <p>Qualitative analysis for: CMV IgG and CMV IgM Hepatitis B surface antigen and core antibody (HBSAG) Hepatitis C (HCV) HIV p24 antigen and antibody to HIV-1 and HIV-2</p> <p>Blood grouping</p> <p>Qualitative analysis for ABO/rhesus factor blood typing Anti-A Anti-B</p>	<p>In house method using Beckman Coulter Access 2 Automated Immunoassay Analyser by chemiluminescence.</p> <p>In house method using Beckman Coulter Access 2 Automated Immunoassay Analyser by chemiluminescence</p> <p>Grifols MD multiscard ABO-D-Rh by Direct Antiglobuin</p>
Human Serum		
Human Serum		
Human whole Blood		



7774
Accredited to
ISO/IEC 17025:2017

Schedule of Accreditation
issued by
United Kingdom Accreditation Service
2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

Bourn Hall Limited, trading as Bourn Hall Clinic
Issue No: 007 Issue date: 17 July 2020

Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
HUMAN BODY FLUIDS Semen	Andrology examination activities for the purposes of clinical diagnosis: <u>Fertility testing</u> Sperm concentration Sperm motility Sperm morphology	Documented in-house procedures based on WHO Laboratory Manual 5 th Edition (2010): Documented in-house procedure Semen Analysis in conjunction with manufacturer's instructions for analysis using CASA Sperminator or manual microscopy for low concentrations Stained slides (KWIK-Diff) microscopy
END		