### **Schedule of Accreditation**

issued by

### **United Kingdom Accreditation Service**

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK



9190

**TQ12 4PB** 

Accredited to ISO/IEC 17025:2017

#### **Axiom Veterinary Laboratories Ltd**

Issue No: 014 Issue date: 11 April 2022

Manor House Contact: Claire Richardson Brunel Road Tel: +44 (0)1626 355 655

Newton Abbot E-Mail: claire.richardson@axiomvetlab.co.uk

Website: http://www.axiomvetlab.com

#### Testing performed at the above address only

#### **DETAIL OF ACCREDITATION**

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
ANIMAL TISSUES and FLUIDS	Biochemistry Tests	
Bovine, Ovine, Serum, Heparinised Plasma	Albumin  ALP AST β Hydroxy butyric acid Total bilirubin Calcium Chloride	Documented In-house Methods CHEM005, CHEM006, CHEM007, CHEM 008 using Olympus AU680
(Serum only) (Serum, OxF plasma)	Creatinine kinase (CK) Creatinine Gamma glutamyl transferase (GGT) Glutamate dehydrogenase (GLDH) Glucose	
(Heparin Blood only)  (Serum only)	Glucose Glutathione peroxidase (GSH-Px) Magnesium Non-esterified fatty acid (NEFA) Inorganic phosphate Potassium Sodium Thyroxine (T4) Total protein Urea	(and in-house method HAEM011)
Ovine Serum, Heparinised Plasma	Vitamin B12	IMM012 using Immulite 2000 Xpi Competitive chemiluminescent enzyme immunoassay

Assessment Manager: KB Page 1 of 4



#### 3130

Accredited to ISO/IEC 17025:2017

### Schedule of Accreditation issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

·

#### **Axiom Veterinary Laboratories Ltd**

Issue No: 014 Issue date: 11 April 2022

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
ANIMAL TISSUES and FLUIDS	Haematology Tests	Documented In-house Methods
Bovine, Ovine, EDTA, Bovine Heparin Blood	Haemoglobin Packed Cell Volume Platelet count Total Red cell count Total White cell count MCH (Mean corpuscular haemoglobin)	Sysmex XT2000i-V Manual HAEM007
	MCV (Mean corpuscular volume)	
	MCHC (Mean corpuscular heamoblobin concentration)	
Ovine, Bovine, EDTA, Blood film	Haemoparasite screen Differential white cell count	HAEM010 & H242 Blood film exam staining & microscopy
ANIMAL TISSUES and FLUIDS	Molecular Tests	Documented In-house Methods
Faeces	Mycobacterium avium subsp. paratuberculosis detection	MDU001-MDU005 using Kingfisher Flex Nucleic Acid Purification System & QuantStudio 5 / ABI 7500 qPCR platform PCR using VetMAX M. paratuberculosis 2.0 kit
Tissue, Serum (single and pooled), Milk (single and bulk)	Detection of Bovine Viral Diarrhoea Virus (BVDV) RNA	MDU001-MDU005 using Kingfisher Flex Nucleic Acid Purification System & QuantStudio 5 / ABI 7500 qPCR platform RT PCR using VetMAX BVDV 4ALL Kit
	Parasitology Tests	
Faeces, Intestinal contents	Detection of Fasciola hepatica (Liver Fluke) and Paramphistomes (Rumen Fluke) eggs	VBA432 Sedimentation

Assessment Manager: KB Page 2 of 4



Accredited to ISO/IEC 17025:2017

#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

#### **Axiom Veterinary Laboratories Ltd**

Issue No: 014 Issue date: 11 April 2022

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used
ANIMAL TISSUES and FLUIDS	Serology Tests	Documented In-house methods using automated DS2 Analyser SEF017 and Dynex Agility SEF025, unless specified
Bovine Serum, Plasma	Detection of antibodies to Bovine Viral Diarrhoea Virus (BVDV)	SEF026 IDEXX BVDV p80 Protein ELISA
Bovine Milk	Detection of antibodies to Bovine Viral Diarrhoea Virus (BVDV)	SEF023 Svanovir BVD-Ab ELISA
Bovine, Ovine, Caprine, Porcine, Equine Serum	Detection of antibodies to Chlamydophila abortus	SEF018 IDScreen Chlamydophila abortus Indirect ELISA
Bovine Milk, Bovine, Ovine, Caprine, Serum, Plasma	Detection of antibodies to Mycobacterium avium subsp. paratuberculosis (Johne's)	SEF032 IDEXX Mycobacterium paratuberculosis (Johnes) ELISA
Bovine, Ovine Serum, Plasma, Milk	Detection of antibodies to Fasciola hepatica (Liver fluke)	SEF010 Bio-XMonoscreen Fasciola Hepatica Ab ELISA
Bovine, Serum	Detection of antibodies to Leptospira Hardjo-bovis	SEF001 PrioCHECK L.hardjo ELISA
Bovine Milk	Detection of antibodies to Leptospira hardjo	SEF022 Linnodee Bovine Leptospira hardjo ELISA
Ovine, Caprine, Serum, Plasma, Milk	Detection of antibodies to Maedivisna Virus (MVV) and caprine arthritis- encephalitis Virus (CAEV)	SEF020 IDVet MVV/CAEV Indirect ELISA
Bovine, Ovine, Caprine, Canine, Equine, Serum, Milk	Detection of antibodies to Neospora canium	SEF013 PrioCHECK Neospora caninum 2.0 Indirect ELISA
Bovine Serum, Plasma	Detection of antibodies to gE Antigen of Infectious Bovine Rhinotracheitis Virus (IBR)	SEF014 ID.VET IBR gE Ab ELISA
Bovine Milk	Detection of antibodies to gE Antigen of Infectious Bovine Rhinotracheitis (IBR)	SEF027 IDEXX IBR (BHV-1) gE ELISA on DS2 only
Bovine Milk	Detection of antibodies to Infectious Bovine Rhinotracheitis (IBR)	SEF024 Svanovir IBR-Ab ELISA
Bovine Serum, Plasma	Detection of antibodies to Infectious Bovine Rhinotracheitis (IBR)	SEF002 IDVet IBR (BOHV-1) IDVet ELISA

Assessment Manager: KB Page 3 of 4



Accredited to ISO/IEC 17025:2017

#### **Schedule of Accreditation** issued by

# United Kingdom Accreditation Service 2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

#### **Axiom Veterinary Laboratories Ltd**

Issue No: 014 Issue date: 11 April 2022

#### Testing performed at main address only

Materials/Products tested	Type of test/Properties measured/Range of measurement	Standard specifications/ Equipment/Techniques used		
ANIMAL TISSUES and FLUIDS (Cont'd)	Serology Tests (Cont'd)	Documented In-house methods using automated DS2 Analyser SEF017 and Dynex Agility SEF025, unless specified (Cont'd)		
Bovine, Ovine, Faeces	Detection of <i>Fasciola hepatica</i> (Liver fluke) Antigen	SEF009 Bio-X Monoscreen F.hepatica Ag ELISA		
Bovine, Ovine, Caprine, Serum, Plasma	Detection of antibodies to Schmallenberg Virus (SBV)	SEF019 IDVet Schmallenberg Virus Competition Multi-species ELISA		
Bovine, Ovine, Caprine, Porcine and Canine, Serum, Plasma	Detection of antibodies to Toxoplasma gondii	SEF021 IDScreen Toxoplasmosis Indirect ELISA Multi-species		
Bovine Serum, Plasma, Whole Blood, Ear notch tissue	Detection of Bovine Viral Diarrhoea Virus (BVDV)	SEF007 IDEXX BVDV Antigen ELISA (Erns ELISA)		
END				

Assessment Manager: KB Page 4 of 4