


# Schedule of Accreditation

issued by

## United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

 <b>26481</b>  Accredited to QSI Standard 2021	<b>Vital Radiology Services Limited</b>  <b>Issue No: 003 Issue date: 25 November 2024</b>	
	<b>Office 3</b> Kylna Business Centre Wood Lane End Hemel Hempstead HP2 4BF	<b>Contact: Sue Daniels</b> <b>Tel: +44 (0)7709634281</b> <b>E-Mail: sue.daniels@vitalrad.co.uk</b> <b>Website: www.vitalrad.co.uk</b>
<b>Examinations performed by the Organisation at the locations specified below</b>		

### Locations covered by the organisation and their relevant activities

#### Service locations:

Primary location details	Activity	Location code
Office 3 Kylna Business Centre Wood Lane End Hemel Hempstead HP2 4BF  <b>Local contact</b> Sue Daniels	Head Office location (No clinical activities)	A

#### Activities performed away from the primary location listed above:

Other location details with full postcode	Activity	Location code
Various locations worldwide (remote working radiologists and radiographers)	TELERADIOLOGY PACS	B



26481

Accredited to  
QSI Standard 2021

## Schedule of Accreditation

issued by

### United Kingdom Accreditation Service

2 Pine Trees, Chertsey Lane, Staines-upon-Thames, TW18 3HR, UK

### Vital Radiology Services Limited

**Issue No:** 003 **Issue date:** 25 November 2024

Examinations/Procedures performed by the Organisation at the locations specified

#### DETAIL OF ACCREDITATION

Modality	Type of Examination/Procedure Performed	Equipment/Techniques used	Location Code
TELERADIOLOGY PACS	Reporting – CT/MRI/CR/Mammography	<p>Teleradiology reporting workstation minimum spec as below:</p> <p><b>Workstation:</b> Intel I5 processor, 16GB Ram, 256GB SSD, Windows 11 OS</p> <p><b>Monitors:</b> Admin monitor (&gt;=21")</p> <p>2 x Diagnostic monitors Eizo Radiforce RX370 or greater– (min 3MP) auto calibrated with RadiNET Pro.</p> <p><b>Dictation device:</b> Philips speech mike II</p> <p><b>PACS Software:</b> Inteleviewer V4-17-1-P395</p>	B
END			